

Marlboro-Chesterfield Pathology, P.C.  
30 Page Street, Pinehurst, NC 28374  
910-687-4188

## Specimen Collection



### Cytolyt Cup

Bladder Washings  
Bronchial Washings  
Cul-de-Sac Fluid  
Cyst Fluid  
Esophageal Brushing  
Kidney Washings  
Pelvic Washings  
Peritoneal Fluid  
Sputum

*Note: For copious amounts of fluid please pour 10cc of fluid into cytolyt and submit the remaining fluid as is.*



### Cytolyt Centrifuge Tube

Fine Needle Aspirations  
Breast  
Cyst Fluid  
Neck Mass Aspirations  
Nipple Discharge  
Parotid  
Thyroid



### Formalin Container

Appendix  
Bladder  
Breast  
Cervical Bx  
Cervix  
Colon Bx  
Colon  
Disc Material  
ECC  
EMBX  
Fallopian Tubes  
Fibroids



### Formalin Container

Gastric Bx  
Hemorrhoids  
Hernia Sac  
Liver  
Placenta  
POC  
Prostate Bx/Chips  
Skin Lesions  
Thyroid  
Toes and Fingers  
Vaginal Tissue  
Vas Deferens



### SurePath Vial

Pap smear Collections



### DNA Collection Device

CT/GC Collection when pap smear is not needed



### 50% Alcohol Soln

Urine Cytology