

Collection Guide

						
	ThinPrep	Aptima	E-Swab	Sterile Cup	Cytolyt	Formalin
Pap	✓					
Trichomonas	✓	✓		✓		
BV Panel	✓	✓				
CV Panel	✓	✓				
CT/NG	✓	✓		✓		
HPV	✓					
HSV I/II	✓	✓				
GBS			✓			
Non-Gynecologic Cytology				✓	✓	
Biopsy						✓
Storage	Ambient	Ambient	Ambient	2-8°C	Ambient	Ambient
Stability	14 Days	14 Days	72 Hours	24 Hours	14 Days	30 Days